

### MODEL T E- 1


physiodispenser offers the implantologist a surgical unit for daily use. The ergonomic design and ease of use ensures fatigue-free operation.

- Ultra-short, powerful electric motor with 60 Ncm torque
- Motor speed range: 50 - 40,000 rpm
- Autosense technology
- Electric motor with cable

## CONTENTS

- CONTROL UNIT ( 230V )
- MOTOR TE 10 WITH 1.8 M CABLE
- Foot control
- Irrigant support
- Universal support
- Mains cable

1 DISPOSABLE IRRIGATING TUBE



  
**autoSense  
TECHNOLOGY**

### MODEL TE-3

**DIGITAL PHYSIODISPENSER** offers the implantologist a surgical unit for daily use. The ergonomic design and ease of use ensures fatigue-free operation.

- Intuitive operating concept, large, easy-to-read display
- 10 Pre Set Programs For Basal And Conventional Implantology
- Ultra-short, powerful electric motor with 60 Ncm torque
- Motor speed range: 50 - 40,000 rpm
- Precise torque control, 5- 60 Ncm
- Automatic thread cutter function
- Electric motor with cable.

## CONTENTS

- CONTROL UNIT ( 230V )
- MOTOR TE 30 WITH 1.8 M CABLE
- Foot control
- Irrigant support
- Mains cable